

TECHNICAL DATA SHEET

## PressPro™ Four Cylinder Small Tonnage Powder Press Machine

*Manual hydraulic press · 15–40 T · pointer or digital gauge · compact manual press*



PressPro™ EPPM4P / EPPM4D series four-cylinder manual powder press — manual hydraulic drive, handwheel and lever, pointer or digital pressure gauge. Image shown for reference; dimensions and configuration vary by model.

---

### CONTENTS

- 01** Product Overview
- 02** Core Advantages
- 03** Specifications
- 04** Applications

### ACS MATERIAL, LLC

959 E Walnut St., Suite 100  
Pasadena, CA 91106, USA

Phone +1 (866) 227-0656

Fax +1 (781) 518-0284

**[contact@acsmaterial.com](mailto:contact@acsmaterial.com)**

Revision Date **09/28/2026**

## SECTION I

## Product Overview

The PressPro™ Four Cylinder Small Tonnage Powder Press Machine (EPPM4P / EPPM4D series) is a compact manual hydraulic laboratory pellet press. The 4 Cylinder Pointer Manual Pellet Press series (EPPM4P) features stylish design, exquisite workmanship and wide applicability. It can be paired with specific molds for use with calcium iron analysis, infrared spectroscopy (IR) and X-ray fluorescence (XRF) analyzers, and is widely applied in research and development fields such as superconductors, batteries, powder metallurgy, advanced ceramics, building materials, cement, and the chemical and pharmaceutical industries. The 4 Cylinder Digital Display Manual Pellet Press series (EPPM4D) is an upgrade of the original series hydraulic press, featuring a digital pressure gauge that directly displays tonnage with a precision of 0.01 t, making it ideal for users requiring high-pressure control accuracy. Four tonnages cover 15–40 T; die molds are optional and not included with the main machine.

## SECTION II

## Core Advantages

- **Two gauge options** — pointer gauge on 4P models, digital display gauge on 4D models showing tonnage directly with 0.01 t precision
- **Manual control, stable holding** — manual pressurization and manual gradual pressurization with pressure stability of  $\leq 1$  MPa / 10 min
- **Four capacity options** — maximum press capacities of 15, 24, 30 or 40 T, each available with a pointer or digital pressure gauge
- **Made for lab pellet work** — presses powders such as KBr, NaCl and Al<sub>2</sub>O<sub>3</sub> into pellets, columns and special-shaped bodies for calcium iron analysis, IR, XRF and materials research

## SECTION III

## Specifications

### Technical data

| Specifications               |                                                                                                           |
|------------------------------|-----------------------------------------------------------------------------------------------------------|
| Brand                        | PressPro™ by ACS Material                                                                                 |
| Models (SKU)                 | EPPM-4P15 / 4D15 · EPPM-4P24 / 4D24 · EPPM-4P30 / 4D30 · EPPM-4P40 / 4D40                                 |
| Series                       | EPPM4P – 4 Cylinder Pointer Manual Pellet Press · EPPM4D – 4 Cylinder Digital Display Manual Pellet Press |
| Available Maximum Capacities | 15 / 24 / 30 / 40 T                                                                                       |
| Pressurization Method        | Manual pressurization / manual gradual pressurization                                                     |
| Pressure Stability           | $\leq 1$ MPa / 10 min                                                                                     |
| Configuration                | Main machine, 1 unit; die molds optional and not included                                                 |

| Specifications                                      |                                                       |                               |                                |                                |
|-----------------------------------------------------|-------------------------------------------------------|-------------------------------|--------------------------------|--------------------------------|
| Parameter                                           | 15 T                                                  | 24 T                          | 30 T                           | 40 T                           |
| Pointer Model                                       | EPPM-4P15                                             | EPPM-4P24                     | EPPM-4P30                      | EPPM-4P40                      |
| Digital Model                                       | EPPM-4D15                                             | EPPM-4D24                     | EPPM-4D30                      | EPPM-4D40                      |
| Pressure Range                                      | 0–15 T / 30 MPa                                       | 0–24 T / 34 MPa               | 0–30 T / 31.5 MPa              | 0–40 T / 30 MPa                |
| Pressure Conversion                                 | 1 MPa = 0.5 T; 1 t = 2 MPa                            | 1 MPa = 0.7 T; 1 t = 1.42 MPa | 1 MPa = 0.95 T; 1 t = 1.05 MPa | 1 MPa = 1.33 T; 1 t = 0.75 MPa |
| Piston Diameter                                     | Ø80 mm                                                | Ø95 mm                        | Ø110 mm                        | Ø130 mm                        |
| Max Piston Stroke                                   | <30 mm                                                | <30 mm                        | <30 mm                         | <50 mm                         |
| Pressure Gauge Display Specification (as published) | 4P: not published · 4D: 0.01 t                        |                               |                                |                                |
| Pressurization Method                               | Manual pressurization / manual gradual pressurization |                               |                                |                                |
| Pressure Stability                                  | ≤1 MPa / 10 min                                       |                               |                                |                                |
| Worktable Diameter                                  | Ø88 mm                                                | Ø108 mm                       | Ø128 mm                        | Ø148 mm                        |
| Workspace                                           | W 80 × 150 mm                                         | W 80 × 150 mm                 | W 92 × 160 mm                  | W 106 × 180 mm                 |
| Overall Dimensions                                  | 245 × 180 × 420 mm                                    | 245 × 180 × 420 mm            | 275 × 210 × 440 mm             | 295 × 230 × 490 mm             |
| Weight                                              | ~42 kg                                                | ~42 kg                        | ~58 kg                         | ~78 kg                         |

Each tonnage is offered as a pointer model (4P) and a digital model (4D) sharing the same press body and specifications. Never exceed the rated pressure of the model concerned. Gauge choice: 4P = pointer gauge; 4D = digital display gauge showing tonnage directly.

#### Scope of delivery

— Standard configuration: EPPM4P / EPPM4D series main machine, 1 unit — Die molds: optional and not included

Final accessories and scope of supply are confirmed in the quotation.

#### SECTION IV

### Applications

Widely used in new materials, superconductors, powder ceramics, new energy sources, batteries, construction materials and metallurgical powders. Compresses various powder materials such as potassium bromide (KBr), sodium chloride (NaCl) and aluminum oxide (Al<sub>2</sub>O<sub>3</sub>); produces pellets, columns and special-shaped bodies for scientific research. Compatible with calcium iron analysis, infrared spectroscopy (IR) and X-ray fluorescence (XRF) analyzers. Suitable for testing the compressive strength of various objects.

Product value: four compact manual tonnages from 15 to 40 T, in pointer or digital display versions — a compact manual hydraulic press for laboratory pellet work.

Revision Date: 09/28/2026

#### DISCLAIMER

PressPro™ is a trademark of ACS Material, LLC. ACS Material, LLC believes that the information in this Technical Data Sheet is accurate and represents the best and most current information available to us. ACS Material makes no representations or warranties, either express or implied, regarding the suitability of the equipment for any application or the accuracy of the information contained within this document. Specifications are typical and may vary; confirm product details and application suitability at the time of quotation. Accordingly, ACS Material will not be responsible for damages resulting from use of or reliance upon this information.