

TECHNICAL DATA SHEET

ALab™ Levo Pipette Controller

Pipette controller · 0.1–100 mL serological pipettes · 3 µm replaceable filter



ALab™ Levo Pipette Controller

The ALab™ Levo pipette controller. Image shown for reference.

CONTENTS

- 01 Product Overview
- 02 Core Advantages
- 03 Specifications
- 04 Applications

ACS MATERIAL, LLC

959 E Walnut St., Suite 100
Pasadena, CA 91106, USA
Phone +1 (866) 227-0656
Fax +1 (781) 518-0284
contact@acsmaterial.com
Revision Date **10/16/2026**

SECTION I

Product Overview

The ALab™ Levo is a pipette controller for glass or plastic serological pipettes (0.1–100 mL) that reduces hands-on time and supports a midway refill.

It absorbs 25 mL of liquid in 4 seconds, weighs 118 g and drains drip by drip with precise control; a 3 µm replaceable hydrophobic filter is fitted.

SECTION II

Core Advantages

- **Fast Suction** — 25 mL in 4 s
- **Precise Drainage** — drip by drip
- **Filter** — 3 µm replaceable
- **Ultra-light** — 118 g
- **Midway Refill** — supported
- **Colors** — five designs

SECTION III

Specifications

Technical data — pipette controller · 3 µm filter · 0.1–100 mL serological pipettes

Specifications	
Brand	ALab™ by ACS Material
Model	Levo
SKU #	ALLEVM00
Product Type	Pipette Controller
Compatible pipettes	Glass or plastic serological pipettes (0.1–100 mL)
Aspiration speed	Can absorb 25 mL liquid in 4 secs
Weight	118 g
Filter	3 µm replaceable hydrophobic filter
Silicone adapter	For 0.1–100 mL pipette tip; fully autoclavable
Colors	Five color designs: blue, green, red, yellow, North Carolina blue

Pipetting steps — Levo pipetting is easy and effortless

Levo pipetting is easy and effortless	
Step	Operation
Manufacturing negative pressure	Squeeze the ball to suck the liquid.
Inhale	The filling rate can be adjusted continuously with the help of the operating rod.
Drainage	Press the operating lever to drain the liquid. The unique valve system ensures that the accuracy of the instrument is not affected by the drainage speed.
Blowout function	When it is necessary to blow out the liquid, when the concave surface reaches the nozzle, press the blow out ball.

Scope of delivery

- **Standard configuration:** one ALab™ Levo pipette controller of the selected color.
- **Additional items:** specify separately in the quotation. The TDS lists the product specifications, not the standard packing list.

SECTION IV

Applications

Suited to pipetting with glass or plastic serological pipettes (0.1–100 mL) that needs fast suction, a midway refill and precise drip-by-drip drainage.

Select the color before ordering.

Revision Date: 10/16/2026

DISCLAIMER

ALab™ is a trademark of ACS Material, LLC. ACS Material, LLC believes that the information in this Technical Data Sheet is accurate and represents the best and most current information available to us. ACS Material makes no representations or warranties, either express or implied, regarding the suitability of the equipment for any application or the accuracy of the information contained within this document. Specifications are typical and may vary; confirm product details and application suitability at the time of quotation. Accordingly, ACS Material will not be responsible for damages resulting from use of or reliance upon this information.