

TECHNICAL DATA SHEET

ALab™ dFlow Electronic Bottle-Top Dispenser

Electronic bottle-top dispenser · 0.1–99.9 mL · 100 µL resolution · 16 speed levels · ISO 8655



ACS Material ALab™ dFlow Electronic Bottle-Top Dispenser

CONTENTS

- 01 Product Overview
- 02 Core Advantages
- 03 Specifications
- 04 Applications

ACS MATERIAL, LLC

959 E Walnut St., Suite 100
Pasadena, CA 91106, USA
Phone +1 (866) 227-0656
Fax +1 (781) 518-0284

contact@acsmaterial.com

Revision Date **08/28/2026**

SECTION I

Product Overview

The ALab™ dFlow Electronic Bottle-Top Dispenser is a motor-driven bottle-top dispenser for reagent delivery in the laboratory. It offers a 0.1–99.9 mL volume range and a 10 mL maximum piston lift as separate specifications, at 100 µL resolution. Dispensing is driven by a stepper motor under digital control. Speed: 16 levels; aspiration and dispensing speed are adjustable. Two operating modes: continuous dispenser and stepper. A TFT display module carries volume selection and control, and the control type is external. A reagent recovery feature reduces waste. Quality Control: ISO 8655. Relative error and coefficient of variation are specified individually at the 10 mL, 5 mL and 1 mL nominal volumes (Section III).

SECTION II

Core Advantages

- Motorized stepper dispensing — smooth, accurate dispensing from a stepper motor under digital motor control
- 16 speed levels — aspiration and dispensing speed are adjustable
- Two operating modes — continuous dispenser and stepper
- Reagent recovery feature — a recovery function to reduce waste

SECTION III

Specifications

Technical data

Specifications	
Volume Range	0.1–99.9 mL
Max. Piston Lift	10 mL
Resolution	100 µL
Volume Accuracy	Specified per nominal volume — see the table below
Speed	16 levels · aspiration and dispensing speed are adjustable
Operating Modes	Continuous dispenser · Stepper
Display / Control Type	TFT display module for volume selection and control · external control
Reagent Recovery	Reagent recovery feature
Operating Temperature Range	15 °C–40 °C
Quality Control	ISO 8655
Media	Bases, saline solutions, acids and some organic solvents
Dimensions · Weight · Power	Confirm at quotation

Volume accuracy — specified per nominal volume

10 mL	$R = 0.5\% \cdot CV = 0.3\%$
5 mL	$R = 1.0\% \cdot CV = 0.6\%$
1 mL	$R = 5.0\% \cdot CV = 3.0\%$

R is relative error and CV is coefficient of variation. Accuracy is specified individually at 10 mL, 5 mL and 1 mL, and the values are not interchangeable: at 1 mL both R and CV are 10 times the corresponding values at 10 mL. Select against the value specified for the volume actually dispensed.

Scope of delivery

- **Sales unit: 1 unit**
- **Control type: external control**
- **Power requirements: confirm at quotation**
- **Package contents and configuration: confirm at quotation**
- **Bottle thread and adapters: confirm at quotation**
- **Options and spare parts: confirm at quotation**

SECTION IV**Applications**

Typical applications for the ALab™ dFlow Electronic Bottle-Top Dispenser include:

- **Bases and saline solutions**
- **Acids — confirm chemical compatibility for the specific reagent and concentration before use**
- **Some organic solvents — suitability depends on the specific reagent and concentration**

Product value: a 0.1–99.9 mL volume range at 100 µL resolution with a 10 mL maximum piston lift, 16 speed levels with adjustable aspiration and dispensing speed, two operating modes, a reagent recovery feature and ISO 8655 quality control, under external control.

Revision Date: 08/28/2026

DISCLAIMER

ALab™ is a trademark of ACS Material, LLC. ACS Material, LLC believes that the information in this Technical Data Sheet is accurate and represents the best and most current information available to us. ACS Material makes no representations or warranties, either express or implied, regarding the suitability of the equipment for any application or the accuracy of the information contained within this document. Specifications are typical and may vary; confirm product details and application suitability at the time of quotation. Accordingly, ACS Material will not be responsible for damages resulting from use of or reliance upon this information.