

TECHNICAL DATA SHEET

ALab™ DispensMate-HF Hydrofluoric Acid Bottle-Top Dispenser

Hydrofluoric acid bottle-top dispenser • 1.0–10 mL • 0.2 mL increment • ceramic piston, cylinder, valve ball and seat • FEP inlet and outlet tubes • ISO 8655



ACS Material ALab™ DispensMate-HF Hydrofluoric Acid Bottle-Top Dispenser

CONTENTS

- 01 Product Overview
- 02 Core Advantages
- 03 Specifications
- 04 Applications

ACS MATERIAL, LLC

959 E Walnut St., Suite 100
Pasadena, CA 91106, USA

Phone +1 (866) 227-0656

Fax +1 (781) 518-0284

contact@acsmaterial.com

Revision Date **08/28/2026**

SECTION I

Product Overview

The ALab™ DispensMate-HF Hydrofluoric Acid Bottle-Top Dispenser is a bottle-top dispenser for hydrofluoric acid, designed for the characteristics of HF. One configuration is listed: a 1.0–10 mL dispensing range with a 0.2 mL increment, complying with ISO 8655. Listed component materials are high-purity ceramic for the piston, piston cylinder, valve ball and valve seat, special alloys for the valve spring, and FEP for the inlet and outlet tubes. Volume is set by a vernier-style adjustable volume guide rail, and a reflux valve design reduces reagent waste and helps prevent leakage during idle periods. The dispenser is designed to be disassembled, cleaned and maintained regularly. Compatible bottle-thread formats listed are S40, GL32, GL38, GL25 and GL28. Hydrofluoric acid is highly corrosive and hazardous to health; handle it only under HF-specific laboratory procedures and the applicable safety data sheet (Section IV).

SECTION II

Core Advantages

- Ceramic piston, cylinder, valve ball and seat — all four in high-purity ceramic materials; the valve spring is a special alloy, not ceramic
- FEP tubing and alloy valve spring — inlet and outlet tubes in FEP, valve spring in special alloys, both selected for corrosion resistance
- Vernier-style volume setting — adjustable volume guide rail for fast, repeatable setting across 1.0–10 mL
- Reflux valve design — reduces reagent waste and helps prevent leakage during idle periods; easy to disassemble and clean

SECTION III

Specifications

Technical data

Specifications	
Type	Bottle-top dispenser for hydrofluoric acid
Dispensing Range	1.0–10 mL · increment 0.2 mL
Working Accuracy	Listed per test volume — see the table below
Standard Compliance	ISO 8655
Piston / Piston Cylinder / Valve Ball / Valve Seat	High-purity ceramic materials
Valve Spring / Tubes	Valve spring in special alloys, selected for corrosion resistance · inlet and outlet tubes in FEP
Volume Setting / Reflux Valve	Vernier-style adjustable volume guide rail · reflux valve reduces reagent waste and helps prevent leakage during idle periods
Compatible Bottle-Thread Formats	S40 · GL32 · GL38 · GL25 · GL28
Dimensions · Weight	Confirm at quotation

Working accuracy — listed per test volume

10 mL	Systematic error 50 µL · 0.5% — random error 20 µL · 0.2%
5 mL	Systematic error 60 µL · 1.2% — random error 20 µL · 0.4%
1 mL	Systematic error 60 µL · 6% — random error 20 µL · 2.0%

Accuracy is listed at three test volumes within the 1.0–10 mL range and the values are not interchangeable: at the 1 mL test volume the listed systematic error is 6% and the random error 2.0%. Accuracy values are provided only for the three test volumes shown. Do not interpolate or extrapolate between the listed points.

Scope of delivery

- ACS sales unit: 1 unit
- Reagent bottle and package configuration: confirm at quotation
- Hydrofluoric acid attacks glass rapidly — confirm bottle and container selection for HF service
- Package contents and configuration: confirm at quotation
- Options and spare parts: confirm at quotation
- Compatible bottle-thread formats are listed in Section III; package configuration is confirmed at quotation

SECTION IV

Applications

Typical applications for the ALab™ DispensMate-HF Hydrofluoric Acid Bottle-Top Dispenser include:

- Dispensing of hydrofluoric acid reagents in trace analysis
- Laboratories where HF is handled under dedicated procedures and containment
-  Hydrofluoric acid is highly corrosive and hazardous to health. Handle it only under your laboratory's HF-specific procedures, with the personal protective equipment, engineering controls and emergency provisions set out in the applicable safety data sheet.
- This dispenser is intended for HF service; it does not remove the need for those controls.

Product value: 1.0–10 mL dispensing in 0.2 mL increments, ISO 8655 standard compliance, listed ceramic, alloy and FEP components, vernier-style volume setting, a reflux valve, and compatibility with the S40, GL32, GL38, GL25 and GL28 bottle-thread formats.

Revision Date: 08/28/2026

DISCLAIMER

ALab™ is a trademark of ACS Material, LLC. ACS Material, LLC believes that the information in this Technical Data Sheet is accurate and represents the best and most current information available to us. ACS Material makes no representations or warranties, either express or implied, regarding the suitability of the equipment for any application or the accuracy of the information contained within this document. Specifications are typical and may vary; confirm product details and application suitability at the time of quotation. Accordingly, ACS Material will not be responsible for damages resulting from use of or reliance upon this information.